

CONSUMABLES

Single and double head high-pressure injectors (83 bar/1,200 psi)

Article number



ELSD 200 ml (Easy Loading Syringe) with filling tube 31 6025-000

Residual volume: 1.5 ml | At 83 bar only usable with pressure sleeve!

ELSD 200 ml (S) (Easy Loading Syringe) 31 6026-000

Residual volume: 1.5 ml | At 83 bar only usable with pressure sleeve!

Packing unit: 50 pieces/carton – 1,500 pieces/pallet



HSD 525 suction tube 31 5080-000

Length: 100 cm
Inner diameter: 2.7 mm
Filling volume: 0.3 ml 3-way stop cock,
6.0 ml tube + spike

Packing unit: 25 pieces/carton

just for
Accutron HP




HSD 527 suction tube with valve 31 5082-000

Length: 100 cm
Inner diameter: 2.7 mm
Filling volume: 0.2 ml valve,
16 ml tube with drip chamber

Packing unit: 50 pieces/carton

just for
Accutron HP



HSD 527 suction tube with valve and inline drip chamber 31 5083-000
(for Scanbag® CM-supply bag)

Length: 100 cm
Inner diameter: 2.7 mm
Filling volume: 0.2 ml valve, 16 ml patient line with drip chamber

Packing unit: 50 pieces/carton

just for
Accutron HP




HDW 526 high-pressure 3-way stop cock 31 5085-100

Filling volume: 0.3 ml

Packing unit: 50 pieces/carton
200 pieces/carton

CONSUMABLES

Single and double head high-pressure injectors (83 bar/1,200 psi)

	Article number
	<p>HS 224 RA patient line</p> <p>Inner diameter: 1.8 mm, Rotating Adapter, Braided PU tube</p> <p>Length: 100 cm, filling volume: 2.5 ml 31 5100-000 Length: 120 cm, filling volume: 3.1 ml 31 5126-000</p> <p>Packing unit: 50 pieces/carton</p>
	<p>HS 224 RA patient line</p> <p>Inner diameter: 1.8 mm, Rotating Adapter, Braided PU tube</p> <p>Length: 150 cm, filling volume: 3.9 ml 31 5136-000 Length: 180 cm, filling volume: 4.6 ml 31 5135-000</p> <p>Packing unit: 25 pieces/carton</p>
	<p>HS 224/180 RA patient line with valve 31 5124-000</p> <p>Inner diameter: 1.8 mm, Rotating Adapter Braided PU tube with check valve</p> <p>Length: 180 cm Filling volume: 4.6 ml</p> <p>Packing unit: 50 pieces/carton</p>
	<p>HS 224 patient line</p> <p>Inner diameter: 1.8 mm Braided PU tube</p> <p>Length: 120 cm, filling volume: 3.1 ml 31 5125-000 Packing unit: 50 pieces/carton</p> <p>Length: 180 cm, filling volume: 4.6 ml 31 5134-000 Packing unit: 25 pieces/carton</p>
	<p>HP patient line</p> <p>Inner diameter: 1.5 mm, PCV-tube</p> <p>Length: 120 cm, Filling volume: 3.8 ml 31 5129-000 Length: 180 cm, Filling volume: 5.7 ml 31 5132-000</p> <p>Packing unit: 50 pieces/carton</p>

Single and double head high-pressure injectors (83 bar/1,200 psi)

Article number



HP RA patient line PVC-flex


Inner diameter: 1.5 mm
Rotating Adapter

Length: 120 cm, Filling volume: 2.18 ml
Packing unit: 50 pieces/carton

31 5130-000

Length: 150 cm, Filling volume: 2.24 ml
Packing unit: 25 pieces/carton

31 5138-000



HP RA patient line with check valve, PVC-flex

Inner diameter: 1.5 mm
Rotating Adapter/PVC-flex tube, with check valve

Length: 120 cm, filling volume: 2.2 ml

31 5137-000

Length: 150 cm, filling volume: 2.7 ml

31 5121-000

Packing unit: 50 pieces/carton



HP 110 RA patient line with check valve, PVC-flex

31 5140-000

Total length: 110 cm
Inner diameter: 1.5 mm
Filling volume: 2.0 ml

Packing unit: 25 pieces/carton

straight packed



Accutron HP-D tube system

31 7183-000

Pressure side	Length	Filling volume
CM-Side	16 cm	0.7 ml
NaCl-Side	16 cm	0.7 ml

Suction Side	Length	Filling volume
CM-Side, green	100 cm	16 ml incl. drip chamber
NaCl-Side, white	100 cm	16 ml incl. drip chamber
2 x drip chamber		each one 10 ml

Packing unit: 25 pieces/carton

just for
Accutron HP-D



ScanBag® HP-D tube system

31 7184-000

(for Scanbag® CM-supply bag)

Pressure side	Length	Filling volume
CM-Side	16 cm	0.7 ml
NaCl-Side	16 cm	0.7 ml

Suction side	Length	Filling volume
CM-Side, green	100 cm	16 ml incl. drip chamber
NaCl-Side, white	100 cm	16 ml incl. drip chamber
2 x drip chamber		each one 10 ml

Packing unit: 25 pieces/carton

just for
Accutron HP-D

CONSUMABLES

Single and double head high-pressure injectors (83 bar/1,200 psi)



Accutron HP-D - Y-Pressure line
with 2 check valves

Article number

31 7187-000

Pressure side	Length	Filling volume
CM-Side	18 cm	0.46 ml
NaCl-Side	18 cm	0.46 ml
Shared line	200 cm	3.5 ml

Packing unit: 25 pieces/carton

just for
Accutron HP-D



Accutron HP-D - Y-connector 83 bar
with 2 CV and RA - long branches

31 7188-000

Pressure side	Length	Filling volume
CM-Side	195 cm	4.85 ml
NaCl-Side	195 cm	4.85 ml
Shared line	15 cm	0.26 ml

Packing unit: 25 pieces/carton

just for
Accutron HP-D



SP 227 vented spike, luer lock

31 9095-100

Filling volume: 0.4 ml

Packing unit: 50 pieces/carton
200 pieces/carton - 11,200 pieces/pallet



Minispike, green

31 9099-100

Filling volume: 0.2 ml

Packing unit: 50 pieces/carton
200 pieces/carton - 11,200 pieces/pallet